

Take a Better Look

5mm and 10mm Laparoscope Line
Bariatric Laparoscope Line
Pediatric Laparoscope Line
Operating Laparoscope Line



Take A Better Look

See Brilliance

See why Stryker is the image leader in high definition imaging with superior contrast and resolution for an unrivaled optical performance. Stryker laparoscopes combine the highest quality sapphires and latest rod lens technology for superior transmission and excellent detail recognition.

Because you demand the best, the improved durability design with laser welded leakage paths and soldered sapphires will withstand repeated autoclave cycles over time.



Bright. Clear. Vivid.

- » High Definition Image
- » Autoclavable
- » Full Screen Picture

Operating Laparoscope Line

- » Working channel for single incision procedures
- » 10mm, 0° Direction of View, 5mm working channel, 45° and 90° ocular angles
- » 11mm, 0° Direction of View, 6mm working channel, 90° ocular angle
- » Compatible with the Filshie® Clip, Falope-Ring® Band, CO₂ lasers

Scope Warmer

Stryker's **Scope Warmer** is a disposable, sterile packaged, self-contained heating device that pre-warms endoscopes to body temperature, overcoming the temperature difference that creates condensation on the endoscope lens.

Stryker's **Scope Warmer** is activated without the need for electrical connections or messy saline solutions.




Scope Warmer
502-360-000 (Box of 6)

Wingman™ Scope Holder

Wingman™ is a pneumatic scope-holding system designed to hold any 5mm or 10mm endoscope and camera.

Extremely effective and easy to use, **Wingman™** volunteers an extra hand in the operating room.



Wingman™ Scope Holder
240-240-999

Ordering Information

5 and 10mm Laparoscope Line

5mm 0° Autoclavable Laparoscope
502-555-010

5mm 30° Autoclavable Laparoscope
502-555-030

5mm 45° Autoclavable Laparoscope
502-555-045

5mm 0° Autoclavable Laparoscope
502-585-010

5mm 30° Autoclavable Laparoscope
502-585-030

5mm 45° Autoclavable Laparoscope
502-585-045

10mm 0° Autoclavable Laparoscope
502-457-010

10mm 30° Autoclavable Laparoscope
502-457-030

10mm 45° Autoclavable Laparoscope
502-457-045

10mm Bariatric Laparoscope Line

10mm 0° Autoclavable Bariatric Laparoscope
502-557-010

10mm 30° Autoclavable Bariatric Laparoscope
502-557-030

10mm 45° Autoclavable Bariatric Laparoscope
502-557-045

Pediatric Laparoscope Line

2.7mm 0° Autoclavable Pediatric Laparoscope
502-543-010

2.7mm 30° Autoclavable Pediatric Laparoscope
502-543-030

Operating Laparoscope Line

11mm, 90° Eyepiece, Operating Laparoscope
331-000-000

10mm, 90° Eyepiece, Operating Laparoscope
331-010-090

10mm, 45° Eyepiece, Operating Laparoscope
331-010-045

Operating Laparoscope Accessories

Detachable Stopcock,
11mm Operating Laparoscope
105-202-754

Detachable Stopcock,
10mm Operating Laparoscope
105-199-070

Port Seal, Small (Blue)
105-189-807

Port Seal, Large (Green)
105-189-808

Literature number: 1000-900-338C
JJ/BL 2.5k 01/07

Copyright © 2007 Stryker
Printed in USA

Operating Laparoscope Accessories (continued)

Laser Coupler
331-110-100

Merrimack Laser Adapter
331-110-110

Coherent Laser Adapter
331-110-120

Sharplan Laser Adapter
331-110-130

Sharplan 5 Laser Adapter
331-110-131

Laparoscope Accessories

10mm Bariatric Laparoscope Sterilization Tray
233-032-111

5mm and 10mm Laparoscope Sterilization Tray
233-032-108

Fog Blocker (Box of 20)
380-022-000

Scope Warmer (Box of 6)
502-360-000

Polishing Paste
105-191-925

Fiberoptic Light Cables

6.5mm x 10ft Fiberoptic Light Cable
(Wolf scope & Stryker light source fittings, snap-on)
233-065-010

5.0mm x 10ft Fiberoptic Light Cable
(Universal scope & light source fittings)
233-050-090

5.0mm x 10ft Clear Case Fiberoptic Light Cable
(Storz scope & Stryker light source fittings, screw-on)
233-050-064

5.0mm x 10ft Fiberoptic Light Cable
(Stryker scope & light source fittings)
233-050-065

5.0mm x 10ft Clear Case Fiberoptic Light Cable
(Wolf scope & Stryker light source fittings, snap-on)
233-050-069

2.0mm x 7.5ft Clear Case Fiberoptic Light Cable
(ACMI scope & Stryker light source fittings)
233-050-067

5.0mm x 10ft Fiberoptic Light Cable
(Storz scope & Stryker light source fittings, screw-on)
233-050-084

Lightpost Adapters

ACMI adapter
331-000-100

Wolf/Dyonics adapter
331-000-200

Storz Adapter
331-000-300

Right Angle adapter
331-000-400

Products referenced with ™ designation are trademarks of Stryker.
Products referenced with ® designation are registered trademarks of Stryker.

stryker®

Joint Replacements

Trauma, Extremities & Deformities

Cranio-maxillofacial

Spine

Biologics

Surgical Products

Neuro & ENT

Interventional Pain

Navigation

Endoscopy

Communications

Imaging

Patient Handling Equipment

EMS Equipment

Rehabilitation Services

5900 Optical Court
San Jose, CA 95138
t: 800.624.4422 f: 800.729.2917

www.stryker.com

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label, and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.